

MADITA-Fun

The modern and integrative therapy chair



MADITA-Fun – the modern and integrative therapy chair

- MADITA-Fun is available in four sizes for children from early support age up to playgroup and school age.
- The possibility of seat height adjustment makes the integration of children in early support, school and at home substantially easier.
- The good adjustability of the seat unit allows secure sitting and provides the child with increased freedom for activities involving the arms.
- The MADITA-Fun therapy chair provides an ideal basis for the construction of individually-moulded seat units.



Seat angle can be adjusted up to 37°



Seat depth adjustment for longer use



The sitting position can easily be lowered via the gas pressure spring from table height to playing height, and allows integrated work with children and siblings of the same age



MADITA-Fun can be equipped with hydraulic height adjustment



Back angle can be adjusted from (-)5° - 25°



The seat width and the abduction can be adjusted through the width, depth and angle adjustable thigh guide



Height and angle adjustable armrests

	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
Back height with short back	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
Back height with long back	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
Seat height	25 - 62 cm	25 - 62 cm	25 - 62 cm	25 - 62 cm
Seat height with electric height adjustment	32 - 62 cm	32 - 62 cm	32 - 62 cm	32 - 62 cm
Seat height with hydraulic height adjustment	27 - 63 cm	27 - 63 cm	27 - 63 cm	27 - 63 cm
Lower leg length	17 - 24 cm	22 - 32 cm	27 - 38 cm	30 - 48 cm
Total width	52 cm	54 cm	64 cm	64 cm
Total length	58 cm	61 cm	75 cm	75 cm
Total width with combination subframe	71 cm	76 cm	86 cm	86 cm
Seat angle	(-)5° - 37°	(-)5° - 37°	(-)5° - 37°	(-)5° - 37°
Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
Max. load	30 kg	30 kg	45 kg	50 kg
Weight	19 kg	20 kg	22 kg	24 kg
HMV No.	26.11.05.0022			

