

MADITA

The mobile therapy chair which grows with the child



MADITA – the mobile therapy chair which grows with the child

- Thanks to the versatile options, sizes and adjustment possibilities, MADITA is for all ages from toddler to teenager.
- Individually adjustable seat dimensions and supporting pelotte pads allow a stable and upright sitting position.
- The seat height can be adjusted to different table or playing heights using a gas pressure spring support.
- Seat tilt is adjustable via gas pressure spring, and can therefore be adapted quickly and easily to the current condition of the child.
- MADITA can be adjusted if required with a slight forward tilt, thus making activities easier and promoting activation of the back musculature.
- During relaxation phases, the therapy chair can be set to horizontal position.
- Special constructions are possible based on the MADITA therapy chair, such as for example the integration of individual body depressions.



Seat tilt is executed with the aid of a gas pressure spring and the height adjustment by foot pedal



Seat and back angle adjustment in the active or passive positions



Seat depth adjustment for longer use



Height and side adjustable thorax pelotte pads



Height and angle adjustable armrests

	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
Back height with short back	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
Back height with long back	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
Seat height	46 - 58 cm	46 - 58 cm	46 - 58 cm	46 - 58 cm
Seat height with low subframe	35 - 40 cm	35 - 40 cm	-	-
Lower leg length	17 - 24 cm	22 - 32 cm	27 - 38 cm	30 - 48 cm
Total width	45 cm	45 cm	53 cm	53 cm
Total length	66 cm	66 cm	77 cm	77 cm
Seat angle	(-)5° - 37°	(-)5° - 37°	(-)5° - 37°	(-)5° - 37°
Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
Max. load	30 kg	30 kg	45 kg	50 kg
Weight	16 kg	17 kg	18 kg	19 kg
HMV No.	26.11.05.0021			

