

MALTE / MALTE Outdoor

The compact and versatile posterior-walkers



MALTE/MALTE Outdoor – the compact and versatile posterior-walkers

- MALTE is a walking aid especially suitable for children from early support up to adolescence.
- MALTE MINI is ultra-light and extremely compact for an early start.
- The abducted frame shape forms a large support surface and at the same time limits pendular motion outwards.
- The individual width and height adjustable handles allow adaptation to the individual body measurements.
- The easy-fold base frame allows problem-free storage and transport in the smallest of spaces (size 0 Mini is not foldable).
- Equipment such as direction locks for the front wheels and individual height and width adjustable handles are already included in the basic model.
- The universally adjustable buttocks pelotte pad allows the correction and support of upright and walking positions and can be turned into the backrest of the seat in one easy step.



The unique folding system makes handling easier



A gauge permits precise height adjustment of the handles



Combined pelvis and buttocks pelotte pad and three-dimensional adjustable forearm supports for stabilisation of the gait pattern



Individually width and height-adjustable handles with brakes



Guide wheel with integrated direction lock



MALTE Outdoor is equipped with a drum brake including a parking brake as standard



The base frame can be folded on sizes 1-4 for easy transport

	MALTE Size 0 MINI*	MALTE Size 1	MALTE Size 2	MALTE Size 3	MALTE Size 4
Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
Height of lower armrests	-	46 - 70 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
Handle spacing	22 - 36 cm	22 - 36 cm	30 - 44 cm	38 - 52 cm	43 - 57 cm
Seat height	-	33 cm	45 cm	60 cm	70 cm
Total width	52 cm	51 cm	61 cm	68 cm	75 cm
Total length	58 cm	58 cm	70 cm	77 cm	87 cm
Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
Max. load	20 kg	30 kg	50 kg	75 kg	100 kg
Weight	3.9 kg	7.1 kg	8.3 kg	8.5 kg	10.2 kg
HMV No.	10.46.02.3029				

* Frame is not foldable

	MALTE Outdoor Size 1	MALTE Outdoor Size 2	MALTE Outdoor Size 3
Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm
Height of lower armrests	46 - 70 cm	63 - 87 cm	79 - 103 cm
Handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm
Seat height	33 cm	45 cm	60 cm
Total width	56 cm	66 cm	73 cm
Total length	66 cm	75 cm	82 cm
Frame height (folded)	25 cm	32 cm	34 cm
Transport dimensions (L x W x H)	66 x 56 x 25 cm	75 x 66 x 32 cm	82 x 73 x 34 cm
Max. load	30 kg	50 kg	70 kg
Weight	7.1 kg	8.3 kg	9.8 kg
HMV No.	10.46.02.3029		