

MARCY

The anterior-walker for more independence and mobility



NEW - now also available in size 3!



MARCY – the anterior-walker for more independence and mobility

- With the coordinated proportions and multitude of options, MARCY allows maximum independence and mobility.
- The handles are individually adjustable in height and width.
- Three-dimensionally adjustable lower armrests are optionally available to stabilise the gait pattern.
- Optional direction locks can be selected for the front wheels for improved directional stability.
- A particular highlight is the option of the adjustable wheelbase, which increases lateral support and permits mobility to be increased step by step using MARCY.
- The size 2+3 base frames can be folded for easy transport.



The wheelbase can be optionally adjusted and set using a gauge




Handles can be singly adjusted in height and width to permit individual adaptation



Three-dimensionally adjustable lower armrests to stabilise the gait pattern



Seat surface



The support surface can be modified by changing the wheelbase



Optional direction locks for the guide wheels



Storage basket



The size 2+3 base frames can be folded for easy transport

	Size 0	Size 1	Size 2	Size 3
 Handle height	30 - 47 cm	38 - 59 cm	53 - 75 cm	63 - 88 cm
Height of lower armrests	41 - 58 cm	49 - 70 cm	63 - 85 cm	73 - 98 cm
Handle spacing	18 - 32 cm	22 - 36 cm	30 - 44 cm	36 - 50 cm
Seat height	25 cm	30 cm	45 cm	45 cm
Total width	44 cm	48 cm	62 cm	71 cm
Total length	36.5 cm	42.5 cm	57.5 cm	61 cm
Total length (telescopic wheelbase)	39 - 47 cm	45 - 55 cm	60 - 74 cm	65 - 79 cm
Max. load	20 kg	30 kg	50 kg	65 kg
Weight	4.2 kg	5 kg	7.4 kg	9.1 kg
HMV No.	10.46.02.3xxx			



Frame colours for size 0+1:  Telemagenta (pink) for size 2:  Lazurite blue metallic for size 3:  Anthracite metallic
 Sky blue  Traffic red  Pure white