

MOMO Therapy bike select

The therapy bike for more balance



MOMO Therapy bike select – the therapy bike for more balance

- The MOMO is a therapy bike which supports in particular the balance training and mobilisation of children and teenagers.
- The suspension force of the stabilizers is easy to adjust and can be adapted to a child's capabilities. The stabilizers can be fixed rigidly to mediate security and to build confidence.
- During the therapy, they can be adjusted to become very gentle suspension stabilizers. The dynamic stabilizers also ensure traction in any riding situation, and balanced travel on bends.
- The diverse range of MOMO series accessories is also compatible with the therapy bike. This permits the provision of customer tailored equipment.
- The frame geometry with the low access point and centre of gravity simplifies handling.
- The diversity of drives and equipment permits individual adaptations.
- The particularly low-set stabilizers permit both dynamic support and rigid adjustment.



The low access point makes independent use and transfer far easier




The selection ranges from a fixed sprocket wheel through to the 7-speed gear hub with reverse option



Therapy bicycle frame with adjustable spring dampers



Therapy bicycle frame as assembly set for standard commercial bicycles (hub installation depth min. 120 mm)



The adjustable suspension intensity and the dynamic stabilizers ensure traction in any riding situation and balanced travel on bends

	12"	16"	20"	24"	26"
Pedal to saddle upper edge*	28 - 35.5 cm	42 - 57.5 cm	48 - 65.5 cm	57 - 74.5 cm	70 - 87.5 cm
Pedal to saddle upper edge**	31 - 39.5 cm	45 - 62 cm	51 - 68 cm	60 - 85 cm	73 - 98 cm
Saddle post tube to front tube	33 cm	39 cm	45 cm	53 cm	58 cm
Wheelbase	64 cm	76 cm	89 cm	104 cm	112 cm
Total length	99 cm	122 cm	148 cm	168 cm	178 cm
Access point height	16 cm	18 cm	20 cm	22 cm	23 cm
Total width	43 cm	52 cm	62 cm	71 cm	76 cm
Max. load	35 kg	50 kg	65 kg	100 kg	100 kg
Dead weight	10.5 kg	13 kg	15 kg	18 kg	20 kg
HMV No.	22.51.01.0xxx				



* with standard saddle post / ** with T-saddle post

- Frame colours for 12": Telemagenta (pink) Sky blue Pure white
- for 16"+20": Melon yellow Marine blue Pure white
- for 24"+26": Anthracite metallic Pure white