

# SOFIE

The individual seat unit that grows with the child



## SOFIE – the individual seat unit that grows with the child

- The comprehensive and diverse range of accessories permits individual equipment and adaptation.
- The seat depth and width as well as the back height can be easily adjusted and are included in the standard functions.
- The seat pad can be selected with anatomical moulding, and thus also allows complex seating arrangements.
- Special constructions are possible on the basis of the SOFIE model, such as for example the integration of individual body depressions.
- SOFIE can be combined with the FLEXI-AKTIV fixation system for support in a sitting position.



Modern and neat design



Back angle adjustment



The SOFIE seat unit can be adapted to all standard subframes (fig. in combination with FREDDY)



Removable back



Seat depth adjustment



The seat width and the abduction can be adjusted through the width, depth and angle adjustable thigh guide



Abducting leg guides

	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
Back height (short back)	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
Back height (long back)	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
Max. load	30 kg	30 kg	45 kg	50 kg
Weight	9 kg	10 kg	10 kg	12 kg
Recommended for subframes	FINN size 0	FINN size 1 + FREDDY size 0	FREDDY size 1 + anti tip bar, FRANKA Size 2, FRIDOLIN 1+2 Size 2	FREDDY size 1 + anti tip bar, FRANKA size 2, FRIDOLIN 1+2 size 2
HMV No.	26.11.04.0xxx			



SOFIE